
Skill Procedures:

Oxygen Administration Chart

I. Usage

1. The following chart shows the various types of oxygen administration devices and their characteristics.

<i>Device</i>	<i>Flow Rate</i>	<i>% of Oxygen Delivered</i>
Nasal Cannula	2-6 liters per minute	24-40%
Pocket Face Mask	10 liters per minute	50-60%
Venturi Mask	Varies by manufacturer	24-50%
Partial Rebreathing Mask	6-10 liters per minute	35-60%
Non-rebreathing Mask	10-15 liters per minute	90%

Skill Procedures:
Oxygen Administration Chart

This Page Intentionally Left Blank
